

## GENETIC CODE WORKSHEET

### Procedure A: Using the first table (identify the protein coded for by the DNA in your flask)

1. Record the flask # \_\_\_\_\_.
2. Record the DNA bases from your flask in ascending numerical order.
3. *Transcribe* the DNA triplets into the corresponding mRNA codons.
4. *Translate* the mRNA codons into amino acids using a Codon Table (from text).
5. Choose your protein from the list (bottom of page) by comparing the 3 amino acids.

DNA triplets	- - -
mRNA codons	- - -
amino acids	(aa1) - (aa2) - (aa3)
name of protein	

### Procedure B: Using the other tables

1. Use one table for each of the remaining proteins. Write the name of the protein on the bottom line of a table then work backwards (from bottom to top of table). Note: there may be more than one mRNA codon and DNA triplet that fits; just choose any one that works.

DNA triplets	- - -
mRNA codons	- - -
amino acids	(aa1) - (aa2) - (aa3)
name of protein	

DNA triplets	- - -
mRNA codons	- - -
amino acids	(aa1) - (aa2) - (aa3)
name of protein	

DNA triplets	- - -
mRNA codons	- - -
amino acids	(aa1) - (aa2) - (aa3)
name of protein	

**insulin****myoglobin****ribonuclease****hemoglobin**

... alanine - serine - valine ...

... arginine - leucine - phenylalanine ...

... serine - glutamine - lysine ...

... leucine - histidine - cysteine ...